4tek S.r.l

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20037 Paderno Dugnano (MI) ITALY

Tel +39 02 91082211 Fax +39 02 91084668

P.IVA C.F. 05133430966

www.4-tek.it

stefano@4-tek.it riccardo@4-tek.it gabriella@4-tek.it

ASPIRATORI AD ANELLO LIQUIDO LIQUID RING ASPIRATION UNIT

MANUALE D'USO / USER MANUAL

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SECTION A – GENERAL FEATURES

- 1 Description**
- 2 Type of Units and selection**
- 3 System's features**
 - 3.1 Principles of operation
 - 3.2 Multiple Units operation
 - 3.3 Installation Criteria
- 4 Installation design**
 - 4.1 Aspiration circuit
 - 4.2 Liquid waste
 - 4.3 Air exhaust
 - 4.4 Water supply
 - 4.5 Electrical wiring
 - 4.6 Amalgam separator

SECTION B – INSTALLATION, USE AND MAINTENANCE

- 5 General Informations**
- 6 Requirements**
 - 6.1 Use
 - 6.2 Noise level
 - 6.3 Ambient temperature
- 7 Installation**
 - 7.1 Transport
 - 7.2 Assembly
 - 7.3 Electrical connection
 - 7.4 Starting
 - 7.5 Technical information for the installer
- 8 Maintenance**
 - 8.1 User's maintenance
 - 8.2 Maintenance by qualified personnel
- 9 Troubleshooting**
- 10 Spare parts and accessories**
- 11 Warranty**

SECTION C – SAFETY INSTRUCTIONS

- 12 Introduction**
 - 12.1 introduction
 - 12.2 Definitions
 - 12.3 Symbols
 - 12.4 Qualified Personnel
- 13 General**
 - 13.1 Basic Rules
 - 13.2 Use
- 14 Installation**
 - 14.1 General safety rules
 - 14.2 Installation
 - 14.3 Electrical connection
 - 14.4 Water and waste lines
 - 14.5 Starting
- 15 Operation**
 - 15.1 General safety rules
- 16 Maintenance and reparation**
 - 16.1 User's maintenance
 - 16.2 General safety rules
 - 16.3 Dismantling
 - 16.4 Reassembling

TABLE 3 - SYMBOLS USED ON THE EQUIPMENT

TABLE 2 – TECHNICAL DATA

SECTION A – GENERAL

1 DESCRIPTION

The liquid ring aspiration system consists of units of which the liquid ring vacuum pump is the essential element.

This type of pump has a cylindrical body in which an eccentric vane impeller rotates inside of it forming a ring of water that gradually increasing and decreasing. This will generate a state of depression, who shall suck the secretions (through cannulas placed in the dental chair) and expel them (through the drain circuit).

During operation the liquid must be continuously fed into the pump, and this is required both to disperse the heat of compression, and to compensate the amount that comes out, among with secretions, through the exhaust.

2 TYPE OF UNITS AND SELECTION

The units described in this manual are the models without cabinet:

OEP11060 – 2 Hp motor pump -1.5 Kw - without control board and without air/water separator, for 4 chairs;

P001 – Like OEP11060 + Control board (KITCNTP001);

P001/S – Like P001 + Air/water separator (S201);

P001/SB – Like P001 / S + Metallic support;

P001/P – Like P001 + Metal silencer (S202);

P001/T – Like P001 + S201 + S202;

P001/TB – Like P001 / T + Metallic support;

OEP10035 – 1 Hp motor pump - 0.75 kW - without control board and without air/water separator, for 2 chairs;

P002 – Like OEP10035 + Control board (KITCNTP002)

P002/S – Like P002 + S201

and the models with cabinet:

UCA1 – Like P001/SB + soundproofing cabinet;

UCA2 – Like P002 / S + soundproofing cabinet;

The model and the number of units can be chosen with the help of next Table 1 and of the par. 3.2., in consideration of the number of chairs to be served at the same time.

WITH 1 UNIT

Chairs working at the same time	Model of the unit	Main pipe		Secondary pipe	
		Diameter mm ³	Max distance ⁴ m	Diameter mm ³	Max distance ⁴ m
2	1 OEP10035 in all its models	32	8	25	6
4	1 OEP11060 in all its models	32	8	32	6

WITH MORE THAN 1 UNITS

Chairs working at the same time	Model of the unit	Main pipe		Secondary pipe	
		Diameter mm ³	Max distance ⁴ m	Diameter mm ³	Max distance ⁴ m
4	2 OEP10035 in all their models	40	20	25	6
8	2 OEP11060 in all their models	40	30	32	6
12	3 OEP11060 in all their models	50	40	32	6
16	4 OEP11060 in all their models	60	50	32	6


Table 1

³ The diameters of the main pipes follow the standard dimensions. The internal dimension of the pipe can have a different of 2-4mm from the dimensions here indicated: this will not cause any difference in the system performance

⁴ With bigger distances, increase the diameter of 10mm for every 5m added

3 SYSTEM'S FEATURE

3.1 Principles of operation

Picture 1

Picture 1 shows the main components of the system and their respective functions.

Once the water supply to the unit is active and power is ON, the unit is put into operation or stopped by the operator, taking or placing the cannula on the tip holder: these actions in fact activates a low tension micro switch (12 V). The water supply is automatic and controlled by a solenoid valve (Art. 7120003C) who opens when the motor starts and closes when the engine stops. This system allows the pump to be activated only when needed, saving water and electricity.

The pump is designed either for heavy loads of work, or for continuous or discontinuous operation. The motor is equipped with a temperature sensor that shuts down operation in the event of overheating due to any cause, including faulty operations.

3.2 Multiple Units operation

When the system is composed of two or more Units, operation can be run in different ways:

A_ *Simple individual*: each Unit serves one or more specific surgeries;

B_ *Simple Alternate*: same as per A, but at the next event the Unit that starts is the one that was idle at the previous event;

C_ *With PLC computer logic*: Units start either as B, but the second Unit starts when more surgeries are in action;

D_ *With vacuum detector* (Art. V003): the first Unit starts when the first surgery is active, the second one starts when the vacuum in the suction line falls below a preset limit; also in this case there is alternation of operation.

3.3 Installation criteria

Units OEP and UCA can be installed in the surgery itself, when there is no better available area, and usually in the cupboard under the sink. In this case it is important be sure that the motor is properly ventilated. Do not place the pump on the floor: it is equipped with special anti-vibration rubbers and its operation is silent and vibration free. In some cases, if you want to absolutely avoid the noise of the Units without cabinet, you can use an additional metallic silencer (Art. S202).

The room in which the Unit is installed must be equipped with:


- a power line adequate to support the load shown on the nameplate of the pump;
- a low-voltage control line from the surgery;
- a water pipe with the right capacity according to the consumption of the pumps;
- a discharge for liquids;
- appropriate dimensions to exhaust air.

In the position deemed most appropriate, between the surgery and the Unit, an amalgam separator (Art. SA650C) can be installed.

4 INSTALLATION DESIGN

4.1 Aspiration Circuit

The suction pipes must be adapted to the diameter of the suction flow. Furthermore, their layout should be done to facilitate the flow of liquids.

Picture 2

The main pipe **P** acts as a collector and has a diameter proportional to the volume of fluid from the surgeries related to it.


The secondary pipes **D**, linking the main line to the various surgeries **R**, have smaller diameters.

Table 1 shows the minimum diameters and maximum lengths of pipes in relation to the number of surgeries and the type of pump used.

Close to the surgeries, at the beginning of pipes, filters should be installed to ensure the flowing in the tubes only to liquid and air.

It is also advisable to install the piping with a slight slope towards the pump so that the system will result in drainage and operation will be benefited.


Suitable inspection traps **B** may be provided at selected locations.

Picture 2.1

4.2 Liquid waste

The liquid discharge pipe must never be raised, but arranged so that the tubes are inclined with a slope of at least 20 mm for 6 m in length. This allows liquid waste to flow into the best possible way, using gravity force. The exhaust circuit of the Unit includes Air/water separator (Art. S201 – Pic 3)

Picture 3

4.3 Air exhaust


The air must be expelled out by a vertical straight pipe with a diameter of 50 mm.

4.4 Water supply

Close to the unit should be provided a proper water outlet to feed the pump, possibly equipped with a shut-off valve, and that will be connected to it via the appropriate tubing and fittings installed on the Unit. In case of multiple units, the necessary fittings should be provided.

4.5 Electrical wiring

The power line, single or three phase, may be adapted to the power of the Unit installed. It's needed consider all the national or local regulations concerning the execution of electrical systems (in particular those relating to circuit grounding and the installation of circuit breakers), as well as for activities involving the work of employees. It is not allowed to use the neutral conductor, even if on the ground, as a protective conductor. In case of multiple units appropriate precautions must be taken.

Picture 4

The control board of the UCA and OEP single phase (one Unit installation) (Pic. 4), is designed to make the user's operating very simple, in fact the starting of the motor and the input of water running in the system are automatically controlled by a low voltage (12 V) micro switch installed near the operator location. For multiple Units or three-phase Units, the low voltage micro switch is at 24 V.

So the electrical wiring must also include a low tension line from the surgery(ies) to the control board(s) of the Unit(s).

4.6 Amalgam separator

The Amalgam separator (Art. SA650C), an optional accessory, can be installed simply hanging it with a wall bracket, and is to be inserted in the aspiration circuit using adequate pipes, usually with a diameter of 32 mm.

SECTION B – INSTALLATION, USE AND MAINTENANCE

5 GENERAL INFORMATIONS

1. All units bear the mark **CE** comply with directive 93/42/CEE MEDICAL DEVICE and following changes (see attached Declaration of Conformity) provided that they are used for the purpose for which they have been manufactured.
2. In order to avoid any danger it is absolutely necessary that all indications in this manual and all other instructions that may supplied including those contained in the Safety Instructions (SECTION C) are adhered to.
3. Custom's tailored models or variants may differ for technical details. Should any problem arise in this connection, please contact 4tek S.r.l.; the model and the Serial Number of the appliance should be given in the enquiry (find them on the plate of the machine) .

6 REQUIREMENTS

6.1 Use

The liquid ring vacuum pumps should be used only with adequate water supply.
Any solids or dirt must be removed before they can enter into the pump, using adequate aspiration filters that are normally included with the tip support, just after the tip.
Avoid, if possible, an outdoor installation; if necessary, contact 4tek Srl for instructions.

6.2 Noise level

Noise levels of the Units are shown in Table 2 - Technical data, on the last page of this manual.

6.3 Ambient temperature

The maximum temperature allowed for the operation of the pumps is +40 ° C.

7 INSTALLATION

7.1 Transport

Make sure that the lifting device has sufficient capacity by referring to the data in Table 2 - Technical data, on the last page of this manual.

7.2 Assembly

To achieve a quiet and vibration-free running it is necessary that a stable base and stable operating conditions be provided.

Make sure the cooling is not prevented and that ventilation openings are not obstructed.

The laying of the pipes must be done safely to avoid any kind of stress on the machine.

The air inlet and outlet may be never obstructed or stopped during operation.

7.3 Electrical connection

Before performing any work on the machine, open the circuit breaker.

Voltage and frequency must match the informations given on the label. Deviations of $\pm 5\%$ of the tension or $\pm 2\%$ of frequency are allowed, because they do not cause a reduction in power.

Connect the ground wire to the proper terminal.

Set the motor thermal protection to match the current value shown on the plate.

Avoid using power converter because high frequency tension and current armonics can cause electromagnetic disturbances and may damage the electronic circuitry.

Nominal current rating for the motor is for an inlet air temperature and an ambient temperature of +40 ° C.

We declare that for tri-phase motors the direction of turning is clockwise, looking the fan's motor from the up side: however the direction is indicated even with an arrow put on the motor.

7.4 Starting

Make the following connections:

Aspiration – Water – Water Discharge – Air Discharge

Connect the power supply to terminals marked "RETE" and the consensus cannula to the appropriate terminals (see Pic 4).

The solenoid valve on the water circuit opens and closes when the pump runs or stops. After being aspirated, air and secretions pass through the filter, reach the pump and then are channelled to the drain. The outlet can be connected to an air / water separator (Art. S201 – Pic 3) equipped with a water recirculator.

Before starting the Unit(s), make sure the shut-off valve for the water supply (see § 4.4) is open. At the first start a few seconds may elapse before the Unit reaches the nominal degree of vacuum, as the whole water circuit has to fill up.

7.5 Technical information for the installer

The units are tested in the test room with the normal water supply network pressure of about 3.5 bar.

If the main pressure is too low, the water system will take some time to be full. During this period starts and stops may be noticed. Leave that happen a few times until the pump will start running regularly.

With the vacuum regulation valve (Art. VRV00000) the desired vacuum generated may be chosen by the operator (usually the Units are supplied with calibration at 40 cmHg)

8 MAINTENANCE

The maintenance required by the unit is minimal. Motor and pump don't require lubrication. It is however recommended to follow the operations described below.

8.1 User's maintenance

- At the end of each work day, clean filters and cannulas with appropriate liquids containing steramine.
- Before stopping the pump at the end of the day, it is recommended to clean the intake and exhaust pipes with chlorine-based cleaners. This simple procedure, which consists in aspirating small amounts of this liquid, can dissolve the calcium deposits in the system and maintain it in perfect hygienic conditions and free of bad smells.
- After surgical operations it is opportune to rinse the machine by sucking in preferably hot water.
- Twice a month check the cleanliness of the containers of the Amalgam Separator (Art. SA650C), where it has been installed as mentioned in § 3.3.
- Once a year check the cleanliness of the water filter.
- External cleaning: the outside of the appliance should be kept free from dust or dirt deposits; pay attention with units without cabinet: they can reach high temperatures.

8.2 Maintenance by qualified personnel *

- Once a year we suggest the replacement of the external pipes, and careful control of the solenoid, discharge and inner tubes.
- If an internal clearing becomes necessary, please ask directly to 4tek Srl, or send the unit to the company for maintenance.

Additional informations, addressed only to qualified personnel, can be obtained directly by 4tek Srl, always showing model and serial number of the pump.

* For the definition of "qualified personnel" see SECTION C - SAFETY in § 12.4

9 TROUBLESHOOTING

Each unit, before leaving the factory, is checked and tested. If, at the beginning or after a certain period of working, problems occur, the following notes can help to identify the causes and possible solutions.

PROBLEM	SOLUTIONS
Motor doesn't run or stops itself	1-2-3-4-5
Motor runs, but vacuum is low even if the vacuum gauge shows right level of vacuum	6-7
Motor runs, but vacuum is low and the vacuum gauge shows wrong level of vacuum	8
Water leaks	9

Possible Solutions:

1. Make sure the circuit breaker (and any other switch) of the surgery is ON.
2. If the unit has a power switch, make sure it is ON.
3. Check if the low voltage switch, near the operator, is ON.
4. If the water supply system is equipped with a check valve, make sure it is open.
5. If the unit is connected directly to the electrical line by a plug, make sure it is inserted into the socket.
6. Where the above actions do not give results, qualified personnel shall inspect all electrical parts.
7. Remove impurities from the filter secretion.
8. Clean filters and suction pipes.
9. Check all lines and pipes starting from the pump to the operator; lock all connections, replace the broken pipes. If the problem persists, contact a technician.
10. If there is water leak, check all connections. If this control does not solve the problem, a technician will check the solenoid to make sure it works properly

10 SPARE PARTS AND ACCESSORIES

All models of 4tek Srl liquid ring vacuum products can be grouped into two sets: one of the models based on the pump OEP11060 (where 110 was in the past the measure of the impeller diameter in mm and 60 mm its height) and the models based on the OEP10035 (where 100 was the measure of the impeller diameter in mm and 35 mm its height).

Each pump is made of 1 bronze kit (body, cap and impeller) and motor with the shaft that fits into the impeller. For OEP11060 the motor is 2 Hp - 1.5 Kw, while for OEP10035 it is 1 Hp - 0.75 kW.

The main spare parts are:

COMPLETE SECRETION FILTER (Art. FLTSECRBR) consisting of:

Body + cover + O-Ring (Art. FLTSEC004C)

Filter cartridge (Art. FLTSEC002)

Their connections.

COMPLETE WATER CIRCUIT (CRA Art [Model]), including:

Water filter (Art. CRA000007)

Pressure regulator (Art. CRA000012)

Pressure gauge 1 / 8 M 0-6 Post Bar D.40 (Art. CRA000008)

Ball valve 1 / 8 F (Art. CRA000013)

Valve (Art. CRA000011)

Solenoid valve G 1 / 8 NC (Art. 7120003C)

Their connections

COMPLETE CONTROL PANEL (Art. KITCNT. ... [Model]), including:

Electric board assembled for P001 and P002 (Art. 7123135B)

Electrical box + wiring

VACUUM REGULATOR VALVE (Art. VRV00000)

VACUUM GAUGE (Art. 7124001) [Only for P001]

CAPACITOR

For P001 36uF (Art. ECO000132)

For P002 25uF (Art. COCS00009/2)

Other spare parts, less common and not on this list, may be provided by 4tek Srl

The main accessories are:

AIR/WATER SEPARATOR (Art. S201)

SILENCER IN IRON (Art. S202)

AMALGAM SEPARATOR (Art. SA650C)

CLAPET VALVE 1" (Art V001) – **CLAPET VALVE 1" ¼** (Art. V001/4)

VACUUM SWITCH (Art. V003)

FROST RESISTANCE HEATING Length 2 m, to twist around the bronze body (Art. RSSCLT2M)

FROST RESISTANCE HEATING DOUBLE Length 4 m, to twist around the bronze body (Art. RSSCLT4M)

11 WARRANTY

All machines and parts made by 4tek Srl are warranted to be free from defects in materials and workmanship for a period of 2 years.

The warranty conditions are described on the card included in unit's documents.

If you need assistance, contact your local representative or directly 4tek Srl

IMPORTANT:

- **damages caused by wrong connection are not warranty case**
- **the unit must be installed and tested by qualified personnel**
- **follow the wiring diagrams provided by the manufacturer to connect the unit**
- **it's important to know the characteristics of electricity supply**
- **it's required to follow the electrical rules with particular reference to what is reported in § 4.5 of this Manual.**

SECTION C - SAFETY INSTRUCTIONS

12 INTRODUCTION

12.1 Introduction

The suction units described in this document are compact and with controlled level of risk. Are however provided the following safety rules that, even if they don't consider all possible cases, must also be followed to meet the regulations.

If particular problems, not fully treated or not treated at all in this manual, occur or if you need more informations, please go directly to 4tek Srl

12.2 Definitions

Installation: includes the design and installation of the unit and its accessories

Exercise: it includes the setting and the use of the unit.

Maintenance: includes control interventions, prevention, investigation of faults and their solutions.


Repairation: includes interventions that involve dismantling, replacement parts and riassembly of the unit.

12.3 Symbols

The following symbols are used with their meaning.

It means that we must pay special attention to technical aspects, which may often not be obvious even to specialists.

To avoid problems that can cause, directly or indirectly, serious personal injuries or property damages, it is equally important to respect all other warnings connected with transport, installation, operation and maintenance.

12.4 Qualified Personnel

Staff who, thanks to training, experience, as well as specific knowledge about relevant rules, is authorized to install, start and make maintenance on the unit, investigating and avoiding accidents.

It is also desirable that this staff has knowledge of first aid, too.

13 GENERAL

13.1 Basic rules:

Because of their electrical and mechanical properties, if the units are not used and maintained in accordance with the requirements, or if subjected to inappropriate interventions, they can be extremely harmful to people or things.

For this reason, it is necessary that the design and execution of mechanical, electrical and idraulic maintenance is implemented and monitored by qualified personnel (with the exception provided in § 16) who needs to know this User Manual.

The correct and safe use of each unit requires appropriate transport and storage, use in the field of intended use and careful maintenance. It must also take into account the informations provided on the unit and on its components.

13.2 Use

The suction unit covered by these instructions are devices used in medical field: so must be used according to their labels and in accordance with any rules listed in the certificates, product documentation and technical documentation. They must also complied with any local laws relevant to the specific area of use.

14 INSTALLATION

14.1 General safety rules

With the exception provided in § 16, safe installation and maintenance of the unit are possible only through specialized personnel.

In particular it's important to follow the requirements regarding the proper use of lifting equipment and tools and use of personal protective equipment (goggles, gloves, etc.).

Keep in mind that after the shutdown of the unit, the moving parts may continue to rotate for several minutes.

14.2 Installation

In bad conditions of working, metal cabinet or pipes can reach high temperatures

It should be considered that air or liquid can be ejected from the outlet or pressure valves at these high temperatures: so it's recommended to position the unit taking care of this. Temperature sensitive elements such as wires or electronic components, should not be near overheat parts of the unit.

Solids and impurities must be removed using suction filters.

The suction unit has cooling and natural ventilation so the air flow must be able to flow inside and outside of it without any stops. The hot air flow should be taken out from the room where the pump is installed.

Be sure that there are no possibility that water can go into the motor.

Avoid any accidental contact with rotating parts.

14.3 Electrical connection

The connecting wires should be adapted to the type of the unit. Connect the unit according to the scheme that is in the box on the control panel or, on machines without control panel, according to the scheme presented in the manual.

This requires that:

- The inside of the box should be clean and free of cables;
- the protective conductor and the earth are connected;
- all screws are tight;
- low minimum distances are respected;
- unused inlets are closed and their caps securely fastened;
- the panel housing is tight
- any metal surface is clean and lightly greased.

14.4 Water and waste line

Before putting in use the unit, connect it to pipes, so any live parts can not be reached by liquid.

All pipes, tanks, valves and fittings must be suitable for material and size to the required pressure and temperature.

According to air discharge pipes, follow the instructions given in this user guide in § 4.

14.5 Starting

If the unit is working without connection to the piping or without fasteners, for example if someone is testing it, it can move suddenly because of the torque of the motor.

For safe operation of the unit it is necessary to respect following conditions:

- installation and operating conditions correspond to the informations on the label and in this manual to: voltage, current, frequency, wiring, mounting, type of protection;
- the unit is properly installed and connected to the pipes;
- the altitude of where you installed the machine has been taken into account when adjusting the pressure regulator valves;
- the possible pipe for the cooling air is not blocked, the cooling effect is not hindered from dirt surface cooling;
- the rotor can rotate without touching any part;
- the machine has the required direction of rotation;
- all fixing screws, connecting elements and electrical conductors are tight;
- earthing and equipotential links in the right way;
- any accessory is connected in the right way and it works properly;
- protective measures were checked against accidental contact with live or moving parts;
- any separate fans is ready start-up and connected in accordance with the direction of rotation required, and, once activated, does not affect the operation of the machine.

The intake and exhaust openings should be arranged so that foreign parts can not be expelled or sucked or blown into the discharge pipe (danger to eyes and skin, risk of poisoning)
When the inlet is open, the inlet should be covered by guards (grids or similar), to prevent aspiration of foreign parts, including parts of the body or clothing.

15 OPERATION

15.1 General safety rules

The covers that prevent contact with live or rotating parts should not be open when the machines are in operation.

The noise of these units is considerably less than 68 dB (A). If it will increase, stop immediately the machine and start the necessary inspections to eliminate the cause of noise.

After prolonged periods of time, you should consider what has been described in § 14.5

16 MAINTENANCE AND REPARATION

In order to identify and eliminate faults before they cause major damages, you should periodically check the unit.

16.1 User's maintenance

Some maintenance can be performed by the user (see § 8) following this manual

16.2 General safety rules

Disconnect the power supply and any accessories before any maintenance or reparation. In addition to the main circuits, you must pay attention to auxiliary circuits (see § 14.2)

Follow this rules:

- Disconnect power supply;
- Do not let it reconnect by itself;
- Check that there is not power supply;

Measures described can be removed only after you have reassembled the machine and end its maintenance

The operating safety of the unit can be guaranteed only if the reparations are carried out using only original spare parts, or at least approved, and if the maintenance instructions are followed correctly

16.3 Dismantling

After removing the screws, parts are no longer held down by centering. Even with a proper removal is therefore not impossible that there is a sudden detachment of heavy parts. This can cause injury and damage. In case of lack of clarity or indecision about what to do, it is recommended to get in touch directly with 4tek Srl

16.4 Reassembling

If requirements imposed by type of IP protection require the use of special sealants, use proper materials, silicon-free and not hardening

TABLE 3 - SYMBOLS USED ON THE EQUIPMENT

	CE marked		Production date
	ON		OFF
	Serial number		Read manual
	Attention: electric shock		

TABLE 2 - TECHNICAL DATA (see last page of this manual)

**TABLE 2 / TABELLA 2
TECHNICAL DATA / DATI TECNICI**

WITHOUT CABINET	Kw	Volt	Hz	r/min	Vac rate: mbar	Capacity L/min	Dim	Kg	Chairs	dB (A)
P001* A richiesta anche 380V	1,5	230	50/60	2800	558,6	750	40X32X42	40	4	66
P001/S	1,5	230	50/60	2800	558,6	750	40X32X42	40	4	66
P001/SB	1,5	230	50/60	2800	558,6	750	55X42X46	45	4	66
P001/P	1,5	230	50/60	2800	558,6	750	40X32X42	40	4	63
P001/T	1,5	230	50/60	2800	558,6	750	40X32X42	40	4	63
P001/TB	1,5	230	50/60	2800	558,6	750	55X42X46	45	4	63
P002	0,75	230	50/60	2800	532	383	32X27X37	35	2	59
P002/S	0,75	230	50/60	2800	532	383	32X27X37	35	2	59
WITH CABINET	Kw	Volt	Hz	r/min	Vac rate: mbar	Capacity L/min	Dim	Kg	Chairs	dB (A)
UCA1	1,5	230	50/60	2800	558,6	750	55X42X46	57	4	61
UCA2	0,75	230	50/60	2800	532	383	40X32X46	45	2	56