CE

# ASPIRATORE AD ANELLO UMIDO

## RAIN1A

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#### **4tek S.r.l**

Viale Dell'Industria 66  
20037 Paderno Dugnano (MI) ITALY  
Tel +39 02 91082211 Fax +39 02 91084668  
P.IVA C.F. 05133430966  
[www.4-tek.it](http://www.4-tek.it)  
[contact@4-tek.it](mailto:contact@4-tek.it)

## Sommario

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## Index

1	WARNING AND CAUTION instructions.....	24
2	Caution.....	25
2.1	To prevent injury .....	25
2.2	Caution : To Reduce Risk of Electric Shock.....	26
2.3	Danger : To Reduce Risk of Explosion or Fire .....	26
2.4	Classification.....	26
3	Installation note .....	27
3.1	Environmental Conditions for normal operation .....	27
3.2	Introduction.....	27
3.3	Unpacking.....	27
3.4	Mounting.....	27
3.5	Installation Location .....	27
3.6	Suction Inlet Hose Connection.....	28
3.7	Draining suction hose connection .....	29
3.8	Exhaust Air Hose Connection .....	29
3.9	Electrical Connection .....	29
3.10	Electrical Shock Hazard .....	29
4	Operation .....	29
4.1	Check List before start up the system.....	29
5	Product Information .....	30
5.1	Functional description.....	30
5.2	Wiring Diagram of Circuit Board .....	32
5.3	Use of Amalgam collecting device .....	34
6	Maintenance.....	35
6.1	Cleaning Suction inlet filter .....	35
6.2	Cleaning Suction Unit .....	35
7	Disposal .....	36
8	Trouble-Shooting Chart .....	37
9	Warranty.....	38

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



No.	Product	Q'ty	Material	Unit	Part No.	Material	Unit	Scale
1	CASE	1	34	LOOK NUT			1	1
2	CHUCK HOOD INCHES	1	38	MESHES		6205Z	1	1
3	CIRCUIT BOARD	1	35	BEARING			1	1
4	POWER SWITCH	1	36	BEARING			1	1
5	CIRCUIT BREAKER	1	37	BRASSING FIBRE PLATE			1	1
6	RUBBER COVER	1	38	O-RING			1	1
7	Cable Guard	2	39	INLET/OUTLET			1	1
8	Cable Guard	1	40	RING BLOWER TOP			1	1
9	RUBBER	4	41	SHIRT			1	1
10	SUCTION INLET (TOP)	1	42	KEY		6b6	1	1
11	O-RING	2	43	BEARING		6202Z	1	1
12	O-RING	2	44	MOTOR UNIT			1	1
13	SUCTION INLET (LOWER)	1	45	MOTOR COVER			1	1
14	O-RING	1	46	COOLING FAN COVER			1	1
15	FILTER	1	47	COOLING FAN COVER			1	1
16	O-RING	1	48	GASKET			4	4
17	SEWAGE OUTLET	1	49	OUTLET COVER			1	1
18	UNION STRAIGHT	1	50	INLET/OUTLET			1	1
19	SEPARATOR	1	51	ELBOW			2	2
20	TOP INLET/OUTLET	1	52	FOOT VALVE			4	4
21	O-RING	1	53	TAPPING BOLT		MAX 12L	2	2
22	TOP INLET/OUTLET CAP	1	54	TAPPING BOLT		MAX 12L	2	2
23	O-RING	1	55	GROSS BOLT		MAX 12BL	10	10
24	LOWER INLET/OUTLET	1	56	HEX GROSS BOLT/NUT		MAX 12BL	4	4
25	SEPARATOR HOUSING TOP	1	57	GROSS BOLT		MAX 12BL	4	4
26	O-RING	1	58	TAPPING BOLT		MAX 12L	6	6
27	SEPARATOR HOUSING LOWER	1	59	SPRING WASHER		W8	1	1
28	AL HOUSING	1	60	SPRING WASHER		W8	1	1
29	OIL SEAL	1	61	TAPPING BOLT		MAX 12BL	3	3
30	O-RING	2	62	GROSS BOLT		MAX 12BL	5	5
31	SEPARATOR PLATE	1	63	HEX SOCKET BOLT/NUT		MAX 12BL	5	5
32	UNION ELBOW	1	64	SWAY GANTRY/ARM AND HOUSING		MAX 12BL	4	4
33	MOTOR	1	65	INLET FIXED SOCKET BOLT		MAX 12BL	1	1
34	HEX SOCKET BOLT	3	66	HEX SOCKET BOLT		MAX 12BL	3	3
35	HEX SOCKET BOLT	3	67	HEX SOCKET BOLT		MAX 12BL	3	3
36	HEX SOCKET BOLT	3	68	HEX SOCKET BOLT		MAX 12BL	3	3
37	HEX SOCKET BOLT	3	69	HEX SOCKET BOLT		MAX 12BL	3	3
38	HEX SOCKET BOLT	3	70	CROSS BOLT		MAX 12BL	4	4
39	CROSS BOLT	4	71	TAPPING BOLT		MAX 12BL	2	2
40	CROSS BOLT	4	72	CROSS BOLT		MAX 12BL	4	4
41	CROSS BOLT	4	73	CROSS BOLT		MAX 12BL	4	4
42	CROSS BOLT	4	74	CROSS BOLT		MAX 12BL	4	4
43	CROSS BOLT	4	75	CROSS BOLT		MAX 12BL	4	4
44	CROSS BOLT	4	76	CROSS BOLT		MAX 12BL	4	4
45	CROSS BOLT	4	77	CROSS BOLT		MAX 12BL	4	4
46	CROSS BOLT	4	78	CROSS BOLT		MAX 12BL	4	4
47	CROSS BOLT	4	79	CROSS BOLT		MAX 12BL	4	4
48	CROSS BOLT	4	80	CROSS BOLT		MAX 12BL	4	4
49	CROSS BOLT	4	81	CROSS BOLT		MAX 12BL	4	4
50	CROSS BOLT	4	82	CROSS BOLT		MAX 12BL	4	4
51	CROSS BOLT	4	83	CROSS BOLT		MAX 12BL	4	4
52	CROSS BOLT	4	84	CROSS BOLT		MAX 12BL	4	4
53	CROSS BOLT	4	85	CROSS BOLT		MAX 12BL	4	4
54	CROSS BOLT	4	86	CROSS BOLT		MAX 12BL	4	4
55	CROSS BOLT	4	87	CROSS BOLT		MAX 12BL	4	4
56	CROSS BOLT	4	88	CROSS BOLT		MAX 12BL	4	4
57	CROSS BOLT	4	89	CROSS BOLT		MAX 12BL	4	4
58	CROSS BOLT	4	90	CROSS BOLT		MAX 12BL	4	4
59	CROSS BOLT	4	91	CROSS BOLT		MAX 12BL	4	4
60	CROSS BOLT	4	92	CROSS BOLT		MAX 12BL	4	4
61	CROSS BOLT	4	93	CROSS BOLT		MAX 12BL	4	4
62	CROSS BOLT	4	94	CROSS BOLT		MAX 12BL	4	4
63	CROSS BOLT	4	95	CROSS BOLT		MAX 12BL	4	4
64	CROSS BOLT	4	96	CROSS BOLT		MAX 12BL	4	4
65	CROSS BOLT	4	97	CROSS BOLT		MAX 12BL	4	4
66	CROSS BOLT	4	98	CROSS BOLT		MAX 12BL	4	4
67	CROSS BOLT	4	99	CROSS BOLT		MAX 12BL	4	4
68	CROSS BOLT	4	100	CROSS BOLT		MAX 12BL	4	4
69	CROSS BOLT	4	101	CROSS BOLT		MAX 12BL	4	4
70	CROSS BOLT	4	102	CROSS BOLT		MAX 12BL	4	4
71	CROSS BOLT	4	103	CROSS BOLT		MAX 12BL	4	4
72	CROSS BOLT	4	104	CROSS BOLT		MAX 12BL	4	4
73	CROSS BOLT	4	105	CROSS BOLT		MAX 12BL	4	4
74	CROSS BOLT	4	106	CROSS BOLT		MAX 12BL	4	4
75	CROSS BOLT	4	107	CROSS BOLT		MAX 12BL	4	4
76	CROSS BOLT	4	108	CROSS BOLT		MAX 12BL	4	4
77	CROSS BOLT	4	109	CROSS BOLT		MAX 12BL	4	4

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
## 1 WARNING AND CAUTION instructions

While reading your manual, please pay close attention to areas labeled :

### WARNING AND CAUTION

**WARNING:** Warnings are given where failure to observe instruction could result in Injury to people.

**CAUTION:** Cautions are found where failure to observe the instruction could result in Damage to the equipment, associated equipment and process

The following symbols of warning will be found on the pump

**Caution / Warning – attention instruction inside of the manual**

**Electric shock – attached label on the unit**

**Hot surface – attached label on the unit**

	Electric shock (Attached label on the product)		Hot surface (Attached label on the product)
	Alternating current		CE mark
	Attention, consulting Accompanying documents		Protective earth (ground)
	"ON" (power: connection from the mains)		"OFF" (power: disconnection from the mains)
	Date of manufacture		Manufacturer
	Serial No		Refer to instruction manual
	Dispose of Electric and Electronic Equipment Correctly		

## 2 Caution

### 2.1 To prevent injury

- Never operate this product if it has a damaged cord or plug. If it is not working properly, has been dropped, damaged or has fallen into water, please contact appropriate service center for examination and repair.
- Keep the cord away from heated surfaces. All electrical products generate heat. To avoid serious burns, NEVER touch suction motor during or immediately after operation. The temperature to suction motor can be reached 120°C.
- Never block any air exhaust valve or place it on a soft surface where the openings may be blocked. The air exhaust valve is for ventilation of the motor inside the housing. Keep the air exhaust valve free of lint, dirt and other foreign objects.

Blocking air exhaust valve can cause the product to overheat and fire

- Never drop or insert fingers or any other object into any openings.
- Use only in well-ventilated areas. The motor on all product are totally enclosed fan cooled.

**WARNING**

Do not operate the pumps in an atmosphere containing flammable or explosive gases/vapors.

- Be sure to properly identify intake and discharge before using machine.
- Protect Unit from contaminants and moisture.

**WARNING**

Do not modify this equipment without authorization of the manufacturer. And excessive use of this equipment may cause malfunction or motor burnt.

### 2.2 *Caution : To Reduce Risk of Electric Shock*

- Do not disassemble. Disassembly or attempted repairs if accomplished Incorrectly can create electrical shock hazard. Refer servicing to qualified service agencies only.
- Do not use this product in or near area where it can fall or be pulled into Water or other liquids.
- Do not reach for this product if it has fallen into liquid. Unplug immediately.
- Never operate this product outdoors in the rain or in a wet area.

### 2.3 *Danger : To Reduce Risk of Explosion or Fire*

- Do not use this unit in or near explosive atmospheres or where aerosol (spray) products are being used.
- Do not use this product near flames.

### 2.4 *Classification*

- Class I according to EU REG 2017/745equipment.
- Ordinary equipment (IPXO).
- Equipment not suitable for use in the presence of a flammable anesthetic Mixture with air or with oxygen or nitrous oxide.
- Continuous operation.

This equipment has been tested and found to comply with the limit for medical devices in Electromagnetic Compatibility EN 55011/A2: 2021, EN 60601-1-2: 2015/A1:2021, IEC61000-3-3: 2013+A1:2017 +A2: 2021

These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity.

However, there is no guarantee that interference will not occur in a particular installation

If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encourage to try to correct

the interference by one or more of the following measures :

- ✓ Reorient or relocate the receiving device.
- ✓ Increase the separation between the equipment.
- ✓ Connect the equipment into an outlet on a circuit different from that to which the other device(s) are connected.
- ✓ Consult that manufacturer or field service technician for help.

### 3 Installation note

**Intended Use:** This product is a dry suction device which generates vacuum pressure by compressing air.

#### 3.1 *Environmental Conditions for normal operation*

- ✓ Temperature : 10°C to 40°C
- ✓ Humidity : 30% to 80% R.H.
- ✓ Atmospheric pressure : 700 to 1060 hPa

#### 3.2 *Introduction*

This manual has been compiled not only for the care and maintenance of 4TEK suction system now in your possession, but as a helpful reference and guide to prevent many problems which can occur if used improperly.

#### 3.3 *Unpacking*

Carefully remove the machine from the shipping case. Preserve all paperwork for future reference. If damage has occurred from shipment a claim must be filed with the carrier immediately; preserve the shipping carton for inspection by the carrier. If you are required to communicate with your dealer or 4TEK, be sure to include your order numbers for quick identification.

#### 3.4 *Mounting*

Rubber feet are attached to the product. Rubber feet are excellent for applications involving a semi-flexible surface such as bench top; they help to isolate noise and eliminate creeping. However, the horizontal position of the motor shaft should be preferred.


#### 3.5 *Installation Location*

The suction system should be located preferably in a clean, dry and well ventilated area. Please be sure not to block the ventilation ports located on the motor housing. The unit should be placed where the surrounding temperature remains between 10°C and 40°C (50°F and 104°F). Always check to insure the location chose is protected from direct or indirect moisture contact.

**Place suction outlet valve higher position than drainage to prevent backflow**

**ATTENTION:** The max distance between the unit and the dental chair has to be 6mt.

## SCHEME CONNECTIONS

Please be sure not to be bended any of suction inlet hose and outlet hose. This can causes losing vacuum, over or damage of suction motor.


Please use exhaust air hose with 40mm diameter hose.

**CAUTION:** Exhaust air temperature increases significantly and discharge drastically too hot for most plastic piping. Therefore, metal piping or heat-resistant pipe is “DANGER-HOT-DO NOT TOUCH.”


### 3.6 Suction Inlet Hose Connection

- Remove the cover of suction inlet and push-in tubing connector in the inlet hole.


### 3.7 Draining suction hose connection

- Use spiral hose internal diameter Ø20mm
- Connect hose to the outlet hole with clockwise direction
- Please be tighten with hose clip or teflon tape
- Always be higher position than drainage to prevent backflow problem

### 3.8 Exhaust Air Hose Connection

- Remove yellow cap from the exhaust air outlet hole.
- Connect exhaust air outlet hose to the hole.
- Widen hose clamp and push in air outlet hose and re-position hose clamp to tighten as follow.

**IMPORTANT:** Always place end of the air outlet hose out of building or window to prevent noisy problem, hot temperature and bad odor.

### 3.9 Electrical Connection

The mains connection must only be made by well-trained electricians. To insure safe installation, check local current and voltage. Serious damage may occur to the motor if it is connected improper voltage.

### 3.10 Electrical Shock Hazard

Disconnect electrical power at the circuit breaker or fuse box before installing this product. Install the system in a location where it will not come in to contact with water or other liquids. Install the system in a location protected from the weather. Electrically ground the system. Failure to follow these instructions can be resulting death, fire or electrical shock.


## 4 Operation

### 4.1 Check List before start up the system

Check if there is any foreign material in the cooling fan part of the motor. From the motor side of the system, verify the blower is rotating in the direction of clockwise. On the system powered by a 3-phase electric power, change the connection of any two (2) wires to reverse blower rotation. Check any abnormal noise or vibration of the suction system. If this happens refer to the trouble shooting guide on this manual. After removing causes of the noise and vibration, re-start the system.

## 5 Product Information

### 5.1 Functional description

1. Cooling fan cover
2. Motor
3. Impeller chamber cover
4. Pump impeller chamber
5. Case
6. Suction inlet hole
7. Filter
8. Draining

9. Amalgam PCB
10. Power On/OFF
11. Overload protective switch
12. PCB
13. Rubber foot
14. Amalgam collecting vessel

## Technical datasheet


### ► Specification

* Main Voltage (VAC)	230
* Frequency (Hz)	50 / 60
* Electric Power (Kw)	0.55 / 0.63
* Vacuum regulated (mbar)	170 / 200
* Suction motor rotation speed (rpm)	2800 / 3440
* No. of treatment units	1unit
* Weight (Kg)	30
* Dimension (L x W x H)	33 x 43 x 66h cm
* Noisy level (dB)	55


## 5.2 Wiring Diagram of Circuit Board

- NOTE
1. CASE (N. 1-9)
  2. SUCTION INLET-SEWAGE OUTLET (N. 10-18)
  3. SEPARATOR (N. 19-32)
  4. MOTOR (N. 33-51)
  5. OTHER PARTS (N. 52-71)
  6. AMALGAM PARTS (N. 72-77)


No.	Product	Standard	Qty	No.	Product	Standard	Qty
32	UNION ELBOW	Ø12	1	64	HEX SOCKET BOLT	MAX10L	1
31	SEPARATOR PLATE	AN241	1	63	HEX SOCKET BOLT	MAX10L	4
30	O-RING	AN241	2	62	CROSS BOLT	MAX10L	5
29	OIL SEAL		1	61	TAPPING BOLT	M3 S20S	3
28	AL HOUSING		1	60	SPRING WASHER	M3	1
27	SEPARATOR HOUSING LOWER		1	59	SPRING WASHER	M3	1
26	O-RING	AN153	1	58	TAPPING BOLT	M3 S12L	6
25	SEPARATOR HOUSING TOP	AN245	1	57	CROSS BOLT	MAX10L	1
24	O-RING	AN245	1	56	HEX CROSS BOLT	MAX10L	4
23	LOWER IMPELLER		1	55	CROSS BOLT	MAX10L	10
22	TOP IMPELLER CAP		1	54	TAPPING BOLT	MAX10L	2
21	O-RING	Ø70	1	53	TAPPING BOLT	MAX12L	4
20	TOP IMPELLER		1	52	FOOT VALVE		2
19	SEPARATOR		1	51	ELBOW		2
18	UNION STRAIGHT	Ø12	1	50	INLET/OUTLET	ISBT	1
17	SEWAGE OUTLET		1	49	OUTLET COVER		1
16	O-RING	S-29	1	48	GASKET		4
15	FILTER		1	47	COOLING FAN COVER		1
14	O-RING	AN239	1	46	COOLING FAN		1
13	SUCTION INLET(LOWER)	S-40	1	45	MOTOR COVER		1
12	O-RING	S-40	2	44	WATER UNIT	Ø20Z	1
11	O-RING	S-44	2	43	BEARING	Ø20Z	1
10	SUCTION INLET (TOP)		1	42	KEY	Ø6	1
9	IMPELLER		4	41	SHAFT		1
8	Cable Grand	PG11	1	40	RING BLOWER TOP		1
7	Cable Grand	PG3	2	39	IMPELLER		1
6	RUBBER COVER		1	38	O-RING		1
5	CIRCUIT BREAKER	7A	1	37	BRANDING PLATE	Ø20Z	1
4	POWER SWITCH		1	36	BEARING	Ø20Z	1
3	CIRCUIT BOARD		1	35	MISHER		1
2	check cover motor		1	34	LOOK UNIT		1
1	CASE	Motomio	1	33	RING BLOWER LOWER		1
	Product	Standard	Qty		Product	Standard	Qty
	Check	Ador	Model		Check	Scale	Model
	ORAM						
	SRMG						
	J. H.						

**ATEK SRL**

RAIN1A  
DM-PL-008

Rev. 0

### 5.3 Use of Amalgam collecting device

#### LED Lamp

- LED lamp is Green on Sensor Normal : Suction is now ready for operation.
- When amalgam collecting vessel is full over 95% yellow lamp is on with beeper sound: press RESET button and change amalgam collecting vessel with new vessel.
- When amalgam collecting vessel is full over 100% red lamp is on with beeper sound: press RESET button and change amalgam collecting vessel with new vessel.

## 6 Maintenance

### **PERIODICAL MAINTENANCE**

Regular maintenance should be performed by dentist or staff according to following terms.

- ✓ Every 6 months to 1 year: checking amalgam collecting vessel and if it is full or dirty, change it.
- ✓ Every 1 year: check the motor status whether it is abnormal. If abnormal bearing sound of the motor is detected or over heating problem happen, call authorized dealer for checking.

Check the plastic hose or joint whether any leakage problem is occur.

If necessary, change all the necessary parts.

- ✓ Every 5 years: replace the motor bearing, rubber part like O-ring, gasket.

#### 6.1 *Cleaning Suction inlet filter*

Cleaning procedures of Filter part

- Pull out filter part from machine in clockwise direction of cover.

- Clean filter in running water for removing particle or debris.

(Recommending it cleaned at least twice per month for avoiding reduction of suction power

- Re-connect filter.

#### 6.2 *Cleaning Suction Unit*

- ✓ After every treatment: For avoiding bad odor and reducing infection possibility, after each patient, suck 1 glass of water before the next treatment.
- ✓ After finishing last treatment of the day: For avoiding bad odor and reducing infection possibility, it is mandatory, also for the purposes of the validity of the guarantee, to use ANTI-FOAMING cleaning & disinfectant bars, in the dental unit filter housing.

#### **CAUTION**

DO NOT use a foaming cleaning & disinfectant agent. It causes backflow of water to the suction motor and can damage it.

## 7 Disposal

### **CAUTION**

This product may be contaminated after using years. Please make the disposal organization aware of this In order that they Can take the proper safety actions.

Uncontaminated plastic components of the product may be recycled.

The built-in electric circuit board and other electric components should be disposed as electric waste.

Other metal components may be disposed of as metallic waste.

The amalgam inside of the vessel should be sent to the authorized company allowed to collect amalgam.

If the product is returned to the local dealer or 4TEK, all connections should be closed so that they are water-tight.

## 8 Trouble-Shooting Chart

This trouble-shooting is carried out by qualified technicians only.

<b>Problem</b>	<b>Reason</b>	<b>Remedy</b>
Unit does not work	Mains supply is not connected Incorrect electrical connection or Power source	Check mains power, fuses and reset Circuit breaker if necessary
Unusual noise increased	Solid particle or debris is coming into the suction motor	Dismantle unit and clean the Impeller of the motor
Excessive vibration	Damaged impeller Motor and/or impeller are dirty	Dismantle unit and clean the Impeller of the motor
Unit is very hot	Wrong wiring Low voltage Suction inlet is clogged Operating pressure or vacuum is too high	Check wiring Supply proper voltage Clean primary filter Install a relief valve and pressure or Vacuum gauge
Suction is too low	Suction inlet is clogged Leaking problem in the suction plumbing	Clean primary filter Check all suction tubes, hoses and connection. Replace if necessary
Suction power is weak & waste water outlet <u>Hole occurring suction</u>	Red rubber drain valve is back- Side position or fitted wrongly	Unscrew 2 bolts & positioned rubber Drain valve correctly. (convex surface Is outside & flat surface is inside)
Motor overload	Low voltage	Check power source Check wire size and wire connections

## 9 Warranty

The devices of 4TEK S.r.l, if installed and sized correctly based on the study, are guaranteed for a period of 24 months or in any case no more than 36 months from the date of issue of the 4TEK invoice.

In order to obtain performance under this warranty, the buyer must promptly (in no event later than 30 days after discovery of the defect) give written notice of the defect to 4TEK SRL or to a representative in your country.

Customer is responsible for freight charges to 4TEK SRL in all cases.

This warranty does not extend to any goods or parts which have been subjected to misuse, lack of maintenance, neglect, damaged by accident or transit damage.

This express warranty excludes all other warranties or representations expressed or implied by any literature, date or person. 4TEK's maximum liability under this exclusive remedy shall never exceed the cost of the subject product. 4TEK reserves the right, at its sole discretion, to refund the purchase price in lieu of repair or replacement.

4TEK is not responsible or liable for indirect or consequential damage of any kind however, including but not limited to those for use of any products, loss of time, inconvenience, lost profit, labor charges, or other incidental or consequential damages with respect to persons, business, or property, whether as a result of breach of warranty, negligence or otherwise.

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



No.	Product	Qty	Material	Weight	Part No.	Standard	Scale
1	CASE	1	34	LOOK NUT			1
2	CIRCUIT BOARD	1	35	MESHES			1
3	CIRCUIT BOARD	1	36	BEARING			1
4	CIRCUIT BOARD	1	38	MESHES			1
5	CIRCUIT BREAKER	1	37	BRASSING FUSED PLATE			1
6	RUBBER COVER	1	38	O-RINGS			1
7	Cable Ground	2	39	INLET/OUTLET			1
8	Cable Ground	1	40	RING BLOWER TOP			1
9	RUBBERS	4	41	SHAFT			1
10	SUCTION INLET (TOP)	1	42	KEY			1
11	O-RING	2	43	BEARING			1
12	O-RING	2	44	MOTOR UNIT			1
13	SUCTION INLET (LOWER)	1	45	MOTOR COVER			1
14	O-RING	1	46	COOLING FAN COVER			1
15	FILTER	1	47	COOLING FAN COVER			1
16	O-RING	1	48	GASKET			4
17	SEWAGE OUTLET	1	49	OUTLET COVER			1
18	UNION STRAIGHT	1	50	INLET/OUTLET			1
19	SEPARATOR	1	51	ELBOW			2
20	TOP INLET/OUTLET	1	52	FOOT VALVE			4
21	O-RING	1	53	TAPPING BOLT			2
22	TOP INLET/OUTLET CAP	1	54	MAX BEL			2
23	O-RING	1	55	CROSS BOLT			10
24	O-RING	1	56	HEX CROSS BOLT/NUT			4
25	SEPARATOR HOUSING TOP	1	57	GROSS BOLT			4
26	O-RING	1	58	TAPPING BOLT			6
27	SEPARATOR HOUSING LOWER	1	59	SPRING WASHER			1
28	AL HOUSING	1	60	SPRING WASHER			1
29	OIL SEAL	1	61	TAPPING BOLT			3
30	O-RING	2	62	GROSS BOLT			5
31	UNION ELBOW	1	63	HEX SOCKET BOLT/NUT			5
32	UNION ELBOW	1	64	SMART GASKET/ASK 40X 40			4
33	MOTOR	1	65	INLET FIXED SOCKET BOLT			1
34	MOTOR	1	66	HEX SOCKET BOLT			3
35	MOTOR	1	67	HEX SOCKET BOLT			4
36	MOTOR	1	68	HEX SOCKET BOLT			3
37	MOTOR	1	69	HEX SOCKET BOLT			4
38	MOTOR	1	70	CROSS BOLT			4
39	MOTOR	1	71	TAPPING BOLT			2
40	MOTOR	1	72	CROSS BOLT			4
41	MOTOR	1	73	CROSS BOLT			4
42	MOTOR	1	74	GATHERING CASE			1
43	MOTOR	1	75	AMALGAM SET			1
44	MOTOR	1	76	AMALGAM SET			1
45	MOTOR	1	77	AMALGAM SET			1
46	MOTOR	1	78	AMALGAM SET			1
47	MOTOR	1	79	AMALGAM SET			1
48	MOTOR	1	80	AMALGAM SET			1
49	MOTOR	1	81	AMALGAM SET			1
50	MOTOR	1	82	AMALGAM SET			1
51	MOTOR	1	83	AMALGAM SET			1
52	MOTOR	1	84	AMALGAM SET			1
53	MOTOR	1	85	AMALGAM SET			1
54	MOTOR	1	86	AMALGAM SET			1
55	MOTOR	1	87	AMALGAM SET			1
56	MOTOR	1	88	AMALGAM SET			1
57	MOTOR	1	89	AMALGAM SET			1
58	MOTOR	1	90	AMALGAM SET			1
59	MOTOR	1	91	AMALGAM SET			1
60	MOTOR	1	92	AMALGAM SET			1
61	MOTOR	1	93	AMALGAM SET			1
62	MOTOR	1	94	AMALGAM SET			1
63	MOTOR	1	95	AMALGAM SET			1
64	MOTOR	1	96	AMALGAM SET			1
65	MOTOR	1	97	AMALGAM SET			1
66	MOTOR	1	98	AMALGAM SET			1
67	MOTOR	1	99	AMALGAM SET			1
68	MOTOR	1	100	AMALGAM SET			1
69	MOTOR	1	101	AMALGAM SET			1
70	MOTOR	1	102	AMALGAM SET			1
71	MOTOR	1	103	AMALGAM SET			1
72	MOTOR	1	104	AMALGAM SET			1
73	MOTOR	1	105	AMALGAM SET			1
74	MOTOR	1	106	AMALGAM SET			1
75	MOTOR	1	107	AMALGAM SET			1
76	MOTOR	1	108	AMALGAM SET			1
77	MOTOR	1	109	AMALGAM SET			1

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Seule la société 4tek S.r.l. peut modifier cette garantie.